

## 1st Molar - Maxillary Single Tubes



Description	Style	Hook	R/L	.018	.022
Standard Edgewise, without distal extension	Weld		R/L		A00210
Standard Edgewise, without distal extension	Bond		R/L		A08012



Standard Edgewise, distal extension	Weld		R/L	A00160	A00101
Standard Edgewise, distal extension	Bond		R/L	A08000	A08010



Standard Edgewise, distal extension long flange	Weld		R/L	A00952	A00953
Standard Edgewise, distal extension long flange	Weld	Mesial	UR/LL	A00954	A00955
Standard Edgewise, distal extension long flange	Weld	Mesial	UL/LR	A00956	A00957
Standard Edgewise, distal extension long flange	Bond	Mesial	UR/LL	A08006	A08016
Standard Edgewise, distal extension long flange	Bond	Mesial	UL/LR	A08007	A08017



FLI® Standard Edgewise	Weld	Mesial	R	A06770	A06774
FLI® Standard Edgewise	Weld	Mesial	L	A06771	A06775
FLI® Standard Edgewise	Bond	Mesial	R	A08770	A08774
FLI® Standard Edgewise	Bond	Mesial	L	A08771	A08775



Long flange, distal offset	Weld		R/L		A00975
Long flange, distal offset	Weld	Mesial	UR/LL		A00977
Long flange, distal offset	Weld	Mesial	UL/LR		A00979



FLI® Tube	Weld	Mesial	R	A06730	A06734
FLI® Tube	Weld	Mesial	L	A06731	A06735
FLI® Tube	Bond	Mesial	R	A08730	A08734
FLI® Tube	Bond	Mesial	L	A08731	A08735







Andrews, convertible	Weld	Mesial	R	A06112	A06152
Andrews, convertible	Weld	Mesial	L	A06113	A06153
Andrews, convertible	Bond	Mesial	R	A08412	A08452
Andrews, convertible	Bond	Mesial	L	A08413	A08453



FLI® One Piece	Bond	Mesial	R	A18750	A18754
FLI® One Piece	Bond	Mesial	L	A18751	A18755








Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	6° D							
0°	0°	6° D							
0°	0°	6° D							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	7° D							
-10°	0°	7° D							
-10°	0°	7° D							
-10°	0°	7° D							
-10°	0°	10° D							
-10°	0°	10° D							

## 1st Molar - Maxillary Single Tubes

Description	Style	Hook	R/L	.018	.022
 Bench Bioprogressive®, convertible Bench Bioprogressive®, convertible Bench Bioprogressive®, convertible Bench Bioprogressive®, convertible	Weld	Mesial	R	<b>A01194</b>	<b>A01294</b>
	Weld	Mesial	L	<b>A01195</b>	<b>A01295</b>
	Bond	Mesial	R	<b>A08294</b>	<b>A08494</b>
	Bond	Mesial	L	<b>A08295</b>	<b>A08495</b>
 Roth, convertible Roth, convertible Roth, convertible Roth, convertible	Weld	Mesial	R	<b>A05664</b>	<b>A05692</b>
	Weld	Mesial	L	<b>A05665</b>	<b>A05693</b>
	Bond	Mesial	R	<b>A08634</b>	<b>A08636</b>
	Bond	Mesial	L	<b>A08635</b>	<b>A08637</b>
 FLI® MBT FLI® MBT FLI® MBT FLI® MBT	Weld	Mesial	R	<b>A05422</b>	<b>A05522</b>
	Weld	Mesial	L	<b>A05423</b>	<b>A05523</b>
	Bond	Mesial	R	<b>A08562</b>	<b>A08566</b>
	Bond	Mesial	L	<b>A08563</b>	<b>A08567</b>
 FLI® Roth FLI® Roth FLI® Roth FLI® Roth	Weld	Mesial	R	<b>A05600</b>	<b>A05602</b>
	Weld	Mesial	L	<b>A05601</b>	<b>A05603</b>
	Bond	Mesial	R	<b>A08200</b>	<b>A08202</b>
	Bond	Mesial	L	<b>A08201</b>	<b>A08203</b>







Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-10°	0°	15° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							

## 1st Molar - Maxillary Double Tubes

Description	Style	Hook	R/L	.018	.022
 Standard Edgewise	Weld		R/L	<b>A00295</b>	
Standard Edgewise	Weld		R/L		<b>A00283</b>
Standard Edgewise	Bond		R/L		<b>A08015</b>
 Standard Edgewise	Weld	M/D	R/L	<b>A00928</b>	
Standard Edgewise	Weld	M/D	R/L		<b>A00929</b>
 Standard Edgewise, distal extension	Weld		R	<b>A00958</b>	<b>A00960</b>
Standard Edgewise, distal extension	Weld		L	<b>A00962</b>	<b>A00964</b>
Standard Edgewise, distal extension	Weld	Mesial	R	<b>A00966</b>	<b>A00968</b>
Standard Edgewise, distal extension	Weld	Mesial	L	<b>A00970</b>	<b>A00972</b>
 Standard Edgewise, convertible	Weld	Mesial	R	<b>A04630</b>	<b>A04632</b>
Standard Edgewise, convertible	Weld	Mesial	L	<b>A04631</b>	<b>A04633</b>
Standard Edgewise, convertible	Bond	Mesial	R	<b>A08030</b>	<b>A08032</b>
Standard Edgewise, convertible	Bond	Mesial	L	<b>A08031</b>	<b>A08033</b>
 Standard Edgewise, distal extension	Weld		R		<b>A00982</b>
Standard Edgewise, distal extension	Weld		L		<b>A00986</b>
Standard Edgewise, distal extension	Weld	Mesial	R		<b>A00990</b>
Standard Edgewise, distal extension	Weld	Mesial	L		<b>A00994</b>
 Wilson®, convertible	Weld	Mesial	R	<b>A04360</b>	<b>A04362</b>
Wilson®, convertible	Weld	Mesial	L	<b>A04261</b>	<b>A04363</b>
Wilson®, convertible	Bond	Mesial	R	<b>A08260</b>	<b>A08262</b>
Wilson®, convertible	Bond	Mesial	L	<b>A08261</b>	<b>A08263</b>
 Combination Maxillary/Mandibular	Weld	Mesial	UL/LR	<b>A06012</b>	
Combination Maxillary/Mandibular	Weld	Mesial	UR/LL	<b>A06013</b>	
Combination Maxillary/Mandibular	Bond	Mesial	UL/LR	<b>A08302</b>	
Combination Maxillary/Mandibular	Bond	Mesial	UR/LL	<b>A08303</b>	
 Combination Maxillary/Mandibular	Weld		R/L	<b>A06014</b>	
Combination Maxillary/Mandibular	Bond		R/L	<b>A08304</b>	

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°				0°	0°	Gingival	.045
0°	0°	0°				0°	0°	Occlusal	.051
0°	0°	0°				0°	0°	Occlusal	.051
0°	0°	0°				0°	0°	Gingival	.045
0°	0°	0°				0°	0°	Occlusal	.051
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	0°				0°	0°	Occlusal	.045
0°	0°	6° D				0°	6° D	Occlusal	.045
0°	0°	6° D				0°	6° D	Occlusal	.045
0°	0°	6° D				0°	6° D	Occlusal	.045
0°	0°	6° D				0°	6° D	Occlusal	.045
0°	0°	6° D				0°	6° D	Gingival	.045
0°	0°	6° D				0°	6° D	Gingival	.045
0°	0°	6° D				0°	6° D	Gingival	.045
0°	0°	6° D				0°	6° D	Gingival	.045
0°	0°	0°				0°	0°	Gingival	.040
0°	0°	0°				0°	0°	Gingival	.040
0°	0°	0°				0°	0°	Gingival	.040
0°	0°	0°				0°	0°	Gingival	.040
0°	0°	0°				0°	0°	Gingival	.040
0°	0°	0°				0°	0°	Gingival	.040

## 1st Molar - Maxillary Double Tubes

Description	Style	Hook	R/L	.018	.022
 Andrews, convertible	Weld	Mesial	R	A06102	A06142
Andrews, convertible	Weld	Mesial	L	A06103	A06143
Andrews, convertible	Bond	Mesial	R	A08402	A08442
Andrews, convertible	Bond	Mesial	L	A08403	A08443
 Andrews, convertible	Weld	Mesial	R	A06106	A06146
Andrews, convertible	Weld	Mesial	L	A06107	A06147
Andrews, convertible	Bond	Mesial	R	A08406	A08446
Andrews, convertible	Bond	Mesial	L	A08407	A08447
 Andrews, convertible not preangled	Weld	Mesial	R	A06104	A06144
Andrews, convertible not preangled	Weld	Mesial	L	A06105	A06145
Andrews, convertible not preangled	Bond	Mesial	R	A08404	A08444
Andrews, convertible not preangled	Bond	Mesial	L	A08405	A08445
 Andrews, convertible not preangled	Weld	Mesial	R	A06100	A06140
Andrews, convertible not preangled	Weld	Mesial	L	A06101	A06141
Andrews, convertible not preangled	Bond	Mesial	R	A08400	A08440
Andrews, convertible not preangled	Bond	Mesial	L	A08401	A08441
 Tip Edge	Weld	Mesial	R		A06080
Tip Edge	Weld	Mesial	L		A06081
 Bench Bioprogressive®, convertible	Weld	Mesial	R	A01178	A01180
Bench Bioprogressive®, convertible	Weld	Mesial	L	A01179	A01181
Bench Bioprogressive®, convertible	Bond	Mesial	R	A08078	A08080
Bench Bioprogressive®, convertible	Bond	Mesial	L	A08079	A08081
 Alexander, convertible	Weld	Mesial	R	A04380	A04390
Alexander, convertible	Weld	Mesial	L	A04381	A04391
Alexander, convertible	Bond	Mesial	R	A08180	A08190
Alexander, convertible	Bond	Mesial	L	A08181	A08191

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-9°	0°	10° D				0°	0°	Gingival	.045
-9°	0°	10° D				0°	0°	Gingival	.045
-9°	0°	10° D				0°	0°	Gingival	.045
-9°	0°	10° D				0°	0°	Gingival	.045
-9°	0°	10° D				0°	0°	Occlusal	.045
-9°	0°	10° D				0°	0°	Occlusal	.045
-9°	0°	10° D				0°	0°	Occlusal	.045
-9°	0°	10° D				0°	0°	Occlusal	.045
-9°	5°	10° D				5°	0°	Occlusal	.045
-9°	5°	10° D				5°	0°	Occlusal	.045
-9°	5°	10° D				5°	0°	Occlusal	.045
-9°	5°	10° D				5°	0°	Occlusal	.045
-9°	5°	10° D				5°	0°	Gingival	.045
-9°	5°	10° D				5°	0°	Gingival	.045
-9°	5°	10° D				5°	0°	Gingival	.045
-9°	5°	10° D				5°	0°	Gingival	.045
-10°	0°	6° D				0°	6° D	Gingival	.036
-10°	0°	6° D				0°	6° D	Gingival	.036
-10°	0°	15° D	0°	0°	0°				
-10°	0°	15° D	0°	0°	0°				
-10°	0°	15° D	0°	0°	0°				
-10°	0°	15° D	0°	0°	0°				
-10°	0°	15° D				0°	0°	Occlusal	.045
-10°	0°	15° D				0°	0°	Occlusal	.045
-10°	0°	15° D				0°	0°	Occlusal	.045
-10°	0°	15° D				0°	0°	Occlusal	.045



## 1st Molar - Maxillary Double Tubes



Description	Style	Hook	R/L	.018	.022
FLI® Roth	Weld	Mesial	R	A05616	A05618
FLI® Roth	Weld	Mesial	L	A05617	A05619
FLI® Roth	Bond	Mesial	R	A08216	A08218
FLI® Roth	Bond	Mesial	L	A08217	A08219
Roth, convertible	Weld	Mesial	R	A05650	A05670
Roth, convertible	Weld	Mesial	L	A05651	A05671
Roth, convertible	Bond	Mesial	R	A08600	A08620
Roth, convertible	Bond	Mesial	L	A08601	A08621
Roth, convertible	Weld	Mesial	R	A05652	A05672
Roth, convertible	Weld	Mesial	L	A05653	A05673
Roth, convertible	Bond	Mesial	R	A08602	A08622
Roth, convertible	Bond	Mesial	L	A08603	A08623
Roth, convertible	Weld	Mesial	R	A05684	A05686
Roth, convertible	Weld	Mesial	L	A05685	A05687
Roth, convertible	Bond	Mesial	R	A08684	A08686
Roth, convertible	Bond	Mesial	L	A08685	A08687
Roth, without cap	Weld	Mesial	R	A05688	A05690
Roth, without cap	Weld	Mesial	L	A05689	A05691
Roth, without cap	Bond	Mesial	R	A08688	A08690
Roth, without cap	Bond	Mesial	L	A08689	A08691

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	-14° D	-14°	0°	14° D				
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Gingival	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D				0° D	14° D	Occlusal	.045
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				
-14°	0°	14° D	-14°	0°	14° D				

## 1st Molar - Maxillary Triple Tubes



Description	Style	Hook	R/L	.018	.022
Ricketts Bioprogressive®, convertible	Weld		R/L	A01120	A01130
Ricketts Bioprogressive®, convertible	Bond		R/L	A08120	A08130
Ricketts Bioprogressive®, convertible	Weld	Mesial	R	A01122	A01132
Ricketts Bioprogressive®, convertible	Weld	Mesial	L	A01123	A01133
Ricketts Bioprogressive®, convertible	Bond	Mesial	R	A08122	A08132
Ricketts Bioprogressive®, convertible	Bond	Mesial	L	A08123	A08133
Ricketts Bioprogressive®	Weld		R/L	A01124	
Ricketts Bioprogressive®	Bond		R/L	A08124	
Ricketts Bioprogressive®	Weld	M/D	R/L	A01125	
Ricketts Bioprogressive®	Bond	M/D	R/L	A08125	
Combination Maxillary/Mandibular	Weld		R/L	A06010	
Combination Maxillary/Mandibular	Bond		R/L	A08300	
Bioprogressive®, convertible	Weld	Mesial	R	A01156	
Bioprogressive®, convertible	Weld	Mesial	L	A01157	
Bioprogressive®, convertible	Bond	Mesial	R	A08056	
Bioprogressive®, convertible	Bond	Mesial	L	A08057	
Bioprogressive®, convertible	Weld	Mesial	R	A04954	
Bioprogressive®, convertible	Weld	Mesial	L	A04955	
Bioprogressive®, convertible	Bond	Mesial	R	A08134	
Bioprogressive®, convertible	Bond	Mesial	L	A08135	
Ricketts Bioprogressive®, convertible	Weld	Mesial	R	A04960	
Ricketts Bioprogressive®, convertible	Weld	Mesial	L	A04961	
Ricketts Bioprogressive®, convertible	Bond	Mesial	R	A08360	
Ricketts Bioprogressive®, convertible	Bond	Mesial	L	A08361	
Bench Trimorphic, convertible	Weld	Mesial	R	A01160	A01170
Bench Trimorphic, convertible	Weld	Mesial	L	A01161	A01171
Bench Trimorphic, convertible	Bond	Mesial	R	A08060	A08070
Bench Trimorphic, convertible	Bond	Mesial	L	A08061	A08071

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	0°	0°	0°	0°	0°	0°	Gingival	.045
0°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
0°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
0°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
0°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
0°	0°	15° D	0°	0°	15° D	0°	0°	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	0°	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	0°	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	0°	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	15° D	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	15° D	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	15° D	Gingival	.045
0°	0°	15° D	0°	0°	15° D	0°	15° D	Gingival	.045
-10°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
-10°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
-10°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045
-10°	0°	15° D	0°	0°	0°	0°	0°	Occlusal	.045

## 1st Molar - Maxillary Triple Tubes



Description	Style	Hook	R/L	.018	.022
SWLF® convertible	Weld	Mesial	R	<b>A05708</b>	<b>A05728</b>
SWLF® convertible	Weld	Mesial	L	<b>A05709</b>	<b>A05729</b>
SWLF® convertible	Bond	Mesial	R	<b>A08532</b>	<b>A08552</b>
SWLF® convertible	Bond	Mesial	L	<b>A08533</b>	<b>A08553</b>



Roth/Trimorphic, convertible	Weld	Mesial	R	<b>A05654</b>	<b>A05674</b>
Roth/Trimorphic, convertible	Weld	Mesial	L	<b>A05655</b>	<b>A05675</b>
Roth/Trimorphic, convertible	Bond	Mesial	R	<b>A08604</b>	<b>A08624</b>
Roth/Trimorphic, convertible	Bond	Mesial	L	<b>A08605</b>	<b>A08625</b>

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-14°	0°	10° D	-14°	0°	10° D	0°	10° D	Occlusal	.045
-14°	0°	10° D	-14°	0°	10° D	0°	10° D	Occlusal	.045
-14°	0°	10° D	-14°	0°	10° D	0°	10° D	Occlusal	.045
-14°	0°	10° D	-14°	0°	10° D	0°	10° D	Occlusal	.045
-14°	0°	14° D	-14°	0°	14° D	0°	14° D	Occlusal	.045
-14°	0°	14° D	-14°	0°	14° D	0°	14° D	Occlusal	.045
-14°	0°	14° D	-14°	0°	14° D	0°	14° D	Occlusal	.045
-14°	0°	14° D	-14°	0°	14° D	0°	14° D	Occlusal	.045









## 1st Molar - Mandibular Single Tubes

	Description	Style	Hook	R/L	.018	.022
	Standard Edgewise, convertible	Weld		R/L	A04840	A04850
	Standard Edgewise, convertible	Bond		R/L	A08040	A08050
	Standard Edgewise, convertible	Weld	Mesial	R	A04842	A04852
	Standard Edgewise, convertible	Weld	Mesial	L	A04843	A04853
	Standard Edgewise, convertible	Bond	Mesial	R	A08042	A08052
	Standard Edgewise, convertible	Bond	Mesial	L	A08043	A08053
 FLI	FLI® Standard Edgewise	Weld	Mesial	R	A06778	A06782
	FLI® Standard Edgewise	Weld	Mesial	L	A06779	A06783
	FLI® Standard Edgewise	Bond	Mesial	R	A08778	A08782
	FLI® Standard Edgewise	Bond	Mesial	L	A08779	A08783
 FLI	FLI® Tube	Weld	Mesial	R	A06766	A06768
	FLI® Tube	Weld	Mesial	L	A06767	A06769
	FLI® Tube	Bond	Mesial	R	A08766	A08768
	FLI® Tube	Bond	Mesial	L	A08767	A08769
 FLI	FLI® MBT	Weld	Mesial	R	A05426	A05526
	FLI® MBT	Weld	Mesial	L	A05427	A05527
	FLI® MBT	Bond	Mesial	R	A08572	A08576
	FLI® MBT	Bond	Mesial	L	A08573	A08577
	SWLF®, convertible	Weld	Mesial	R	A05714	A05734
	SWLF®, convertible	Weld	Mesial	L	A05715	A05735
	SWLF®, convertible	Bond	Mesial	R	A08538	A08558
	SWLF®, convertible	Bond	Mesial	L	A08539	A08559
	Alexander, convertible not preangled	Weld	Mesial	R	A04384	A04394
	Alexander, convertible not preangled	Weld	Mesial	L	A04385	A04395
	Alexander, convertible not preangled	Bond	Mesial	R	A08184	A08194
	Alexander, convertible not preangled	Bond	Mesial	L	A08185	A08195
	Ricketts Bioprogressive®, without cap	Weld	Mesial	R	A04968	
	Ricketts Bioprogressive®, without cap	Weld	Mesial	L	A04969	
	Ricketts Bioprogressive®, without cap	Bond	Mesial	R	A08368	
	Ricketts Bioprogressive®, without cap	Bond	Mesial	L	A08369	

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-20°	0°	0°							
-20°	0°	0°							
-20°	0°	0°							
-20°	0°	0°							
-20°	0°	8° D							
-20°	0°	8° D							
-20°	0°	8° D							
-20°	0°	8° D							
-22°	-6°	5° D							
-22°	-6°	5° D							
-22°	-6°	5° D							
-22°	-6°	5° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							



## 1st Molar - Mandibular Single Tubes

Description	Style	Hook	R/L	.018	.022
 Ricketts Bioprogressive®, convertible	Weld	Mesial	R	<b>A04978</b>	
Ricketts Bioprogressive®, convertible	Weld	Mesial	L	<b>A04979</b>	
Ricketts Bioprogressive®, convertible	Bond	Mesial	R	<b>A08378</b>	
Ricketts Bioprogressive®, convertible	Bond	Mesial	L	<b>A08379</b>	
 FLI® Tube	Weld	Mesial	R	<b>A06740</b>	<b>A06744</b>
FLI® Tube	Weld	Mesial	L	<b>A06741</b>	<b>A06745</b>
FLI® Tube	Bond	Mesial	R	<b>A08740</b>	<b>A08744</b>
FLI® Tube	Bond	Mesial	L	<b>A08741</b>	<b>A08745</b>
 FLI® One Piece	Bond	Mesial	R	<b>A18758</b>	<b>A18762</b>
FLI® One Piece	Bond	Mesial	L	<b>A18759</b>	<b>A18763</b>
 Bench Bioprogressive®, without cap	Weld	Mesial	R	<b>A01196</b>	
Bench Bioprogressive®, without cap	Weld	Mesial	L	<b>A01197</b>	
Bench Bioprogressive®, without cap	Bond	Mesial	R	<b>A08296</b>	
Bench Bioprogressive®, without cap	Bond	Mesial	L	<b>A08297</b>	
 Bench Bioprogressive®, convertible	Weld	Mesial	R	<b>A01198</b>	
Bench Bioprogressive®, convertible	Weld	Mesial	L	<b>A01199</b>	
Bench Bioprogressive®, convertible	Bond	Mesial	R	<b>A08298</b>	
Bench Bioprogressive®, convertible	Bond	Mesial	L	<b>A08299</b>	
 Roth, convertible, not preangled	Weld	Mesial	R	<b>A05666</b>	<b>A05694</b>
Roth, convertible, not preangled	Weld	Mesial	L	<b>A05667</b>	<b>A05695</b>
Roth, convertible, not preangled	Bond	Mesial	R	<b>A08638</b>	<b>A08640</b>
Roth, convertible, not preangled	Bond	Mesial	L	<b>A08639</b>	<b>A08641</b>
 Andrews, convertible, not preangled	Weld	Mesial	R	<b>A06120</b>	<b>A06160</b>
Andrews, convertible, not preangled	Weld	Mesial	L	<b>A06121</b>	<b>A06161</b>
Andrews, convertible, not preangled	Bond	Mesial	R	<b>A08420</b>	<b>A08460</b>
Andrews, convertible, not preangled	Bond	Mesial	L	<b>A08421</b>	<b>A08461</b>
 FLI Roth, not preangled	Weld	Mesial	R	<b>A05608</b>	<b>A05610</b>
FLI Roth, not preangled	Weld	Mesial	L	<b>A05609</b>	<b>A05611</b>
FLI Roth, not preangled	Bond	Mesial	R	<b>A08208</b>	<b>A08210</b>
FLI Roth, not preangled	Bond	Mesial	L	<b>A08209</b>	<b>A08211</b>






Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-24°	-5°	12° D							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	5° D							
-25°	0°	5° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-30°	-1°	0°							
-30°	-1°	0°							
-30°	-1°	0°							
-30°	-1°	0°							
-30°	+2°	0°							
-30°	+2°	0°							
-30°	+2°	0°							
-30°	+2°	0°							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							

## 1st Molar - Mandibular Double Tubes

Description	Style	Hook	R/L	.018	.022
 Bioprogressive®, Maxillary/Mandibular	Bond	Center	R/L	A08020	
 Wilson®	Weld	Mesial	R	A04370	A04372
 Wilson®	Weld	Mesial	L	A04371	A04373
 Standard Edgewise, convertible	Weld		R/L	A04836	A04846
 Standard Edgewise, convertible	Bond		R/L	A08036	A08046
 Standard Edgewise, convertible	Weld	Mesial	R	A04838	A04848
 Standard Edgewise, convertible	Weld	Mesial	L	A04839	A04849
 Standard Edgewise, convertible	Bond	Mesial	R	A08038	A08048
 Standard Edgewise, convertible	Bond	Mesial	L	A08039	A08049
 Bioprogressive®, not preangulated	Weld	Center	R/L	A01190	
 Bioprogressive®, not preangulated	Bond	Center	R	A08392	
 Bioprogressive®, not preangulated	Bond	Center	L	A08393	
 Ricketts Bioprogressive	Weld	Mesial	R	A01032	
 Ricketts Bioprogressive	Weld	Mesial	L	A01033	
 Ricketts Bioprogressive	Bond	Mesial	R	A08114	
Ricketts Bioprogressive	Bond	Mesial	L	A08115	
Bioprogressive®, convertible	Weld	Mesial	R	A04934	
Bioprogressive®, convertible	Weld	Mesial	L	A04935	
Bioprogressive®, convertible	Bond	Mesial	R	A08034	
Bioprogressive®, convertible	Bond	Mesial	L	A08035	
Ricketts Bioprogressive®, convertible	Weld	Mesial	R	A04964	
Ricketts Bioprogressive®, convertible	Weld	Mesial	L	A04965	
Ricketts Bioprogressive®, convertible	Bond	Mesial	R	A08364	
Ricketts Bioprogressive®, convertible	Bond	Mesial	L	A08365	
Ricketts Bioprogressive®, convertible	Weld	Mesial	R	A04986	
Ricketts Bioprogressive®, convertible	Weld	Mesial	L	A04987	
Ricketts Bioprogressive®, convertible	Bond	Mesial	R	A08386	
Ricketts Bioprogressive®, convertible	Bond	Mesial	L	A08387	

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	0°	0°	0°	0°	0°				
0°	-5°	0°	0°	-5°	0°				
0°	-5°	0°	0°	-5°	0°				
0°	-5°	0°	0°	-5°	0°				
-22°	-5°	12° D	-22°	-5°	12° D				
-22°	-5°	12° D	-22°	-5°	12° D				
-22°	-5°	12° D	-22°	-5°	12° D				
-22°	-5°	12° D	-22°	-5°	12° D				
-22°	-5°	12° D	0°	0°	0°				
-22°	-5°	12° D	0°	0°	0°				
-22°	-5°	12° D	0°	0°	0°				
-22°	-5°	12° D	0°	0°	0°				
-24°	-5°	12° D	-24°	-5°	12° D				
-24°	-5°	12° D	-24°	-5°	12° D				
-24°	-5°	12° D	-24°	-5°	12° D				
-24°	-5°	12° D	-24°	-5°	12° D				
-24°	-5°	12° D	0°	0°	0°				
-24°	-5°	12° D	0°	0°	0°				
-24°	-5°	12° D	0°	0°	0°				
-24°	-5°	12° D	0°	0°	0°				

## 1st Molar - Mandibular Double Tubes

Description	Style	Hook	R/L	.018	.022
 SWLF®, convertible	Weld	Mesial	R	A05712	A05732
SWLF®, convertible	Weld	Mesial	L	A05713	A05733
SWLF®, convertible	Bond	Mesial	R	A08536	A08556
SWLF®, convertible	Bond	Mesial	L	A08537	A08557
 Tip Edge	Weld	Mesial	R		A06082
Tip Edge	Weld	Mesial	L		A06083
 Bench Bioprogressive®, convertible	Weld	Mesial	R	A01164	A01174
Bench Bioprogressive®, convertible	Weld	Mesial	L	A01165	A01175
Bench Bioprogressive®, convertible	Bond	Mesial	R	A08064	A08074
Bench Bioprogressive®, convertible	Bond	Mesial	L	A08065	A08075
 FLI® Roth, not preangled	Weld	Mesial	R	A05620	A05622
FLI® Roth, not preangled	Weld	Mesial	L	A05621	A05623
FLI® Roth, not preangled	Bond	Mesial	R	A08220	A08222
FLI® Roth, not preangled	Bond	Mesial	L	A08221	A08223
 Roth/Andrews, convertible not preangled	Weld	Mesial	R	A05658	A05678
Roth/Andrews, convertible not preangled	Weld	Mesial	L	A05659	A05679
Roth/Andrews, convertible not preangled	Bond	Mesial	R	A08608	A08628
Roth/Andrews, convertible not preangled	Bond	Mesial	L	A08609	A08629










## 1st Molar - Mandibular Triple Tubes

Description	Style	Hook	R/L	.018	.022
 Roth/Andrews, not preangled	Weld	Mesial	R	A05660	A05680
Roth/Andrews, not preangled	Weld	Mesial	L	A05661	A05681
Roth/Andrews, not preangled	Bond	Mesial	R	A08610	A08630
Roth/Andrews, not preangled	Bond	Mesial	L	A08611	A08631

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-25°	0°	8° D	-25°	0°	8° D				
-25°	0°	8° D	-25°	0°	8° D				
-25°	0°	8° D	-25°	0°	8° D				
-25°	0°	8° D	-25°	0°	8° D				
-27°	0°	0°				0°	0°	Gingival	.036
-27°	0°	0°				0°	0°	Gingival	.036
-27°	-5°	6° D	0°	0°	0°				
-27°	-5°	6° D	0°	0°	0°				
-27°	-5°	6° D	0°	0°	0°				
-27°	-5°	6° D	0°	0°	0°				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				
-30°	-1°	4° D	-30°	-1°	4° D				

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-30°	-1°	4° D	-30°	-1°	4° D	-1°	4° D	Occlusal	.045
-30°	-1°	4° D	-30°	-1°	4° D	-1°	4° D	Occlusal	.045
-30°	-1°	4° D	-30°	-1°	4° D	-1°	4° D	Occlusal	.045
-30°	-1°	4° D	-30°	-1°	4° D	-1°	4° D	Occlusal	.045







## 2nd Molar - Maxillary Single Tubes

Description	Style	Hook	R/L	.018	.022
 FLI® Standard Edgewise	Weld	Mesial	R	<b>A06772</b>	<b>A06776</b>
FLI® Standard Edgewise	Weld	Mesial	L	<b>A06773</b>	<b>A06777</b>
FLI® Standard Edgewise	Bond	Mesial	R	<b>A08772</b>	<b>A08776</b>
FLI® Standard Edgewise	Bond	Mesial	L	<b>A08773</b>	<b>A08777</b>
 FLI® Mini Tube	Bond	Mesial	R	<b>A08800</b>	<b>A08806</b>
FLI® Mini Tube	Bond	Mesial	L	<b>A08801</b>	<b>A08807</b>
Bioprogressive®	Weld	Mesial	R	<b>A04982</b>	
Bioprogressive®	Weld	Mesial	L	<b>A04983</b>	
 Bioprogressive®	Bond	Mesial	R	<b>A08382</b>	
Bioprogressive®	Bond	Mesial	L	<b>A08383</b>	
Andrews	Weld	Mesial	R	<b>A06110</b>	<b>A06150</b>
Andrews	Weld	Mesial	L	<b>A06111</b>	<b>A06151</b>
 Andrews	Bond	Mesial	R	<b>A08410</b>	<b>A08450</b>
Andrews	Bond	Mesial	L	<b>A08411</b>	<b>A08451</b>
Andrews not preangulated	Weld	Mesial	R	<b>A06108</b>	<b>A06148</b>
Andrews not preangulated	Weld	Mesial	L	<b>A06109</b>	<b>A06149</b>
 Andrews not preangulated	Bond	Mesial	R	<b>A08408</b>	<b>A08448</b>
Andrews not preangulated	Bond	Mesial	L	<b>A08409</b>	<b>A08449</b>
FLI®	Weld	Mesial	R	<b>A06732</b>	<b>A06736</b>
FLI®	Weld	Mesial	L	<b>A06733</b>	<b>A06737</b>
 FLI®	Bond	Mesial	R	<b>A08732</b>	<b>A08736</b>
FLI®	Bond	Mesial	L	<b>A08733</b>	<b>A08737</b>
 FLI® Mini Tube	Bond	Mesial	R	<b>A08802</b>	<b>A08808</b>
FLI® Mini Tube	Bond	Mesial	L	<b>A08803</b>	<b>A08809</b>
Alexander	Weld	Mesial	R	<b>A04382</b>	<b>A04392</b>
Alexander	Weld	Mesial	L	<b>A04383</b>	<b>A04393</b>
 Alexander	Bond	Mesial	R	<b>A08182</b>	<b>A08192</b>
Alexander	Bond	Mesial	L	<b>A08183</b>	<b>A08193</b>
 FLI® One Piece	Bond	Mesial	R	<b>A18752</b>	<b>A18756</b>
FLI® One Piece	Bond	Mesial	L	<b>A18753</b>	<b>A18757</b>

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	12° D							
0°	0°	12° D							
0°	0°	12° D							
0°	0°	12° D							
-9°	0°	10° D							
-9°	0°	10° D							
-9°	0°	10° D							
-9°	0°	10° D							
-9°	5°	10° D							
-9°	5°	10° D							
-9°	5°	10° D							
-9°	5°	10° D							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	6° D							
-10°	0°	6° D							
-10°	0°	6° D							
-10°	0°	6° D							
-10°	0°	10° D							
-10°	0°	10° D							












## 2nd Molar - Maxillary Single Tubes

Description	Style	Hook	R/L	.018	.022
 Bench Bioprogressive® Trimorphic	Weld	Mesial	R	A01162	A01172
Bench Bioprogressive® Trimorphic	Weld	Mesial	L	A01163	A01173
Bench Bioprogressive® Trimorphic	Bond	Mesial	R	A08062	A08072
Bench Bioprogressive® Trimorphic	Bond	Mesial	L	A08063	A08073
 FLI® Mini Tube	Bond	Mesial	R	A08804	A08810
FLI® Mini Tube	Bond	Mesial	L	A08805	A08811
 FLI® MBT	Weld	Mesial	R	A05424	A05524
FLI® MBT	Weld	Mesial	L	A05425	A05525
FLI® MBT	Bond	Mesial	R	A08564	A08568
FLI® MBT	Bond	Mesial	L	A08565	A08569
 SWLF®, convertible	Weld	Mesial	R	A05710	A05730
SWLF®, convertible	Weld	Mesial	L	A05711	A05731
SWLF®, convertible	Bond	Mesial	R	A08534	A08554
SWLF®, convertible	Bond	Mesial	L	A08535	A08555
 FLI® Roth	Weld	Mesial	R	A05604	A05606
FLI® Roth	Weld	Mesial	L	A05605	A05607
FLI® Roth	Bond	Mesial	R	A08204	A08206
FLI® Roth	Bond	Mesial	L	A08205	A08207
 Roth Trimorphic	Weld	Mesial	R	A05656	A05676
Roth Trimorphic	Weld	Mesial	L	A05657	A05677
Roth Trimorphic	Bond	Mesial	R	A08606	A08626
Roth Trimorphic	Bond	Mesial	L	A08607	A08627

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-10°	0°	12° D							
-10°	0°	12° D							
-10°	0°	12° D							
-10°	0°	12° D							
-14°	0°	0°							
-14°	0°	0°							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	10° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							
-14°	0°	14° D							

## 2nd Molar - Mandibular Single Tubes

Description	Style	Hook	R/L	.018	.022
 FLI® Standard Edgewise	Weld	Mesial	R	<b>A06780</b>	<b>A06784</b>
FLI® Standard Edgewise	Weld	Mesial	L	<b>A06781</b>	<b>A06785</b>
FLI® Standard Edgewise	Bond	Mesial	R	<b>A08780</b>	<b>A08784</b>
FLI® Standard Edgewise	Bond	Mesial	L	<b>A08781</b>	<b>A08785</b>
 FLI® Mini Tube	Bond	Mesial	R	<b>A08812</b>	<b>A08818</b>
FLI® Mini Tube	Bond	Mesial	L	<b>A08813</b>	<b>A08819</b>
 FLI® MBT	Weld	Mesial	R	<b>A05428</b>	<b>A05528</b>
FLI® MBT	Weld	Mesial	L	<b>A05429</b>	<b>A05529</b>
FLI® MBT	Bond	Mesial	R	<b>A08574</b>	<b>A08578</b>
FLI® MBT	Bond	Mesial	L	<b>A08575</b>	<b>A08579</b>
 FLI® Mini Tube	Bond	Mesial	R	<b>A08814</b>	<b>A08820</b>
FLI® Mini Tube	Bond	Mesial	L	<b>A08815</b>	<b>A08821</b>
 FLI® Tube	Weld	Mesial	R	<b>A06742</b>	<b>A06746</b>
FLI® Tube	Weld	Mesial	L	<b>A06743</b>	<b>A06747</b>
FLI® Tube	Bond	Mesial	R	<b>A08742</b>	<b>A08746</b>
FLI® Tube	Bond	Mesial	L	<b>A08743</b>	<b>A08747</b>
 FLI® Mini Tube	Bond	Mesial	R	<b>A08816</b>	<b>A08822</b>
FLI® Mini Tube	Bond	Mesial	L	<b>A08817</b>	<b>A08823</b>
 FLI® One Piece	Bond	Mesial	R	<b>A18760</b>	<b>A18764</b>
FLI® One Piece	Bond	Mesial	L	<b>A18761</b>	<b>A18765</b>
 Bench Bioprogressive® Trimorphic	Weld	Mesial	R	<b>A01166</b>	<b>A01176</b>
Bench Bioprogressive® Trimorphic	Weld	Mesial	L	<b>A01167</b>	<b>A01177</b>
Bench Bioprogressive® Trimorphic	Bond	Mesial	R	<b>A08066</b>	<b>A08076</b>
Bench Bioprogressive® Trimorphic	Bond	Mesial	L	<b>A08067</b>	<b>A08077</b>
 Alexander, not preangled	Weld	Mesial	R	<b>A04386</b>	<b>A04396</b>
Alexander, not preangled	Weld	Mesial	L	<b>A04387</b>	<b>A04397</b>
Alexander, not preangled	Bond	Mesial	R	<b>A08186</b>	<b>A08196</b>
Alexander, not preangled	Bond	Mesial	L	<b>A08187</b>	<b>A08197</b>

Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
0°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-10°	0°	0°							
-14°	0°	0°							
-14°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	0°							
-25°	0°	5° D							
-25°	0°	5° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-5°	6° D							
-27°	-6°	6° D							
-27°	-6°	6° D							
-27°	-6°	6° D							
-27°	-6°	6° D							

## 2nd Molar - Mandibular Single Tubes



Description	Style	Hook	R/L	.018	.022
FLI® Roth, not preangled	Weld	Mesial	R	<b>A05612</b>	<b>A05614</b>
FLI® Roth, not preangled	Weld	Mesial	L	<b>A05613</b>	<b>A05615</b>
FLI® Roth, not preangled	Bond	Mesial	R	<b>A08212</b>	<b>A08214</b>
FLI® Roth, not preangled	Bond	Mesial	L	<b>A08213</b>	<b>A08215</b>
Roth, not preangled	Weld	Mesial	R	<b>A05662</b>	<b>A05682</b>
Roth, not preangled	Weld	Mesial	L	<b>A05663</b>	<b>A05683</b>
Roth, not preangled	Bond	Mesial	R	<b>A08612</b>	<b>A08632</b>
Roth, not preangled	Bond	Mesial	L	<b>A08613</b>	<b>A08633</b>
Ricketts Bioprogressive®	Weld	Mesial	R	<b>A04990</b>	
Ricketts Bioprogressive®	Weld	Mesial	L	<b>A04991</b>	
Ricketts Bioprogressive®	Bond	Mesial	R	<b>A08390</b>	
Ricketts Bioprogressive®	Bond	Mesial	L	<b>A08391</b>	
Andrews, not preangulated	Weld	Mesial	R	<b>A06126</b>	<b>A06166</b>
Andrews, not preangulated	Weld	Mesial	L	<b>A06127</b>	<b>A06167</b>
Andrews, not preangulated	Bond	Mesial	R	<b>A08426</b>	<b>A08466</b>
Andrews, not preangulated	Bond	Mesial	L	<b>A08427</b>	<b>A08467</b>



Main Arch Slot			Auxiliary Slot			Headgear Slot			
Torque	Angulation	Offset	Torque	Angulation	Offset	Angulation	Offset	Location	Size
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-30°	-1°	4° D							
-34°	-5°	6° D							
-34°	-5°	6° D							
-34°	-5°	6° D							
-34°	-5°	6° D							
-35°	+2°	0°							
-35°	+2°	0°							
-35°	+2°	0°							
-35°	+2°	0°							